## Characterization of acrylics for partial and full dentures and Individualization of confectioned acrylic teeth

Finish and polish the denture with paper discs. Sand blast the teeth with aluminium oxide (50 micron). Thoroughly clean the denture with an ultrasound cleaner or a steam cleaner.

Ensure that the sections to be painted are not touched after cleaning. The surface to be painted is moistened with Primer liquid and cured in a light curing unit.

**Note!** If the surface is non-sticky and shiny, the curing time was too long. If there are moist spots of Primer left, the curing time was too short. The colors are put onto a mixing palette. After the characterization of the denture, it must be cured in a light curing unit. An intermediate polymerization of each step is unnecessary. The dispersion layer may be removed with a soft cloth soaked with a bit of Stains-Cleaner. The concluding high-polish is done conventionally with polishing paste and linen buff.

### Primer (thinner) use

The primer could be used to dilute colours in case more fluid consistency is requested. The Primer has to be used before the Stains application, when the upper layers are to be applied, in comparison of small thin lines of staining, in case stains are not to be used inside the material, but outside. When used outside it is suggested to abrade the surface of the restoration using the aluminium oxide sand of 50 micron before the primer application.

The Primer is light cured by the UVA light (200-260 nm or 320-340 nm)(for example using Spektra 2000 for 1 min., using Spektramat for 5 min.)

### Storage

- Do not store with temperature over 25°C. Use the product on the ambient temperature. Shake before using. Close well the top of the container and always close after use. Avoid the direct sunlight exposure.
- Do not use the product after the expiration date indicated on the label.
- Medical device, for dental use only: keep away from children.

### Notes

- Using instructions were realized in base of our experiences and they are to be considered as suggestions for the consumer.
- The producer follows continue research and could bring material variants to improve the product.
- The product has to be used as suggested. The producer doesn't have any responsibility for the damages caused by incorrect use.
- When used for subjects with high sensibility, could cause a product sensibilization. Don't use this product in case of allergic subjects or subjects with well known reactions for some of the components (for example Methylmethacrylate).
- Avoid the direct contact with not cured material. Check all the informations specified on the security sheets (available on request).

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## Ena Tender Stain

### Light curing stains

#### Light curing stains on methacrylate base for dentists and dental technicians

They are available in the following colours:

white, blue, brown, brown2, black, yellow, orange, khaki, red.

- 1 thinner (primer)
- 1 cleaner

#### Main Components

Organic matrix: Methylmethacrylate, Bis-GMA

Additives: densifier agents, catalysts, stabilizers, pigments

#### Indications

Characterization of composite and acrylic materials.

#### Contraindications and collateral effects

This product has been specially developed for the above mentioned indications and is only to be processed as described above. The user is obliged to check before use, whether the product is adequate for use in the intended manner. A sensitization by the product may not be ruled out completely. The use of the product should be discontinued, if any reactions should occur or if there is a known intolerance against one or several of the components (e. g. methylmetacrylate).

#### Instructions for Stains use

Be careful during the light curing stains application because of material fluidity. It's possible to use stains mixed with composite bodies for shoulder and dentin. You can obtain the individualization of the composite mixing, applying or introducing colours. When handling the material be careful not to create bubbles.

If more fluid consistency is required, the colour is to be mixed with Thinner. The brush C is to be used for stains application. The M brush is suggested for composite application and smoothing, while the F brush was formed for the fossa molding.

PROCESSING PHASES: according to the light conditions ca. 5 min.

#### CURING PHASES:

- Laboratory high powerful light curing units  
(as Spektra 2000 - Spektramat) ca. 3 - 5 min.
- Laboratory 36W light curing units  
(as LAMPADA PLUS with light 71 all colours  
as LAMPADA with light 78 only light colours) ca. 30 min.
- Dentist study light curing units  
(according to: kind of unit, colours, thickness) ca. 1-3 min.

Note: For curing you must use a light curing appliance with protected spectrum among 350 nm and 500 nm. It's possible to reach a satisfactory hardening only by using the light curing units with not damaged walls or stoves. Moreover it is necessary to provide for a regular control of the light following producers instructions.