ENA CEM HF

(GB) ENGLISH **pag. 2**

EnaCem HF instructions

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EnaCem HF istruzioni

(D) DEUTSCH **pag.12**

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(F) FRANCAIS **pag.17**

EnaCem HF instructions

(E) ESPAÑOL **pag.22**

EnaCem HF instrucciones

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ENA CEM HF

(GB) ENA CEM HF

EnaCem HF is a dual curing radiopaque fluorescent luting composite, in dentine colour (UD1,UD2,UD3,UD4), for cementation of posts and core build-ups, for luting of ceramic and composite inlays, onlays veneers and crowns. It follows the standard ISO 4049:2000.

Same colours are available in a Try-in paste for colour testing in luting procedures.

Composition

- MONOMER MATRIX: Diurethandimethacrylate; 1,4 - Butandioldimethacrylate.
- TOTAL CONTENT OF FILLERS: 77 % weight; glass filler: mean particle size 4,3 µm and 0.7 µm; highly dispersed silicone dioxide: mean particle size 0,04 µm

Indications. EnaCem HF is used for :

- Cementation of Post
- Cementation of Ceramic laminated veneers, inlays, onlays, jacket crowns, crowns and bridges
- Cementation of Composite veneers, inlays, onlays, crowns and bridges
- Other adhesive lutings
- Core build-up
- Liner for posterior direct and indirect composite restorations

Contra-indications

Uncured resin could cause skin allergy: User should use gloves. In case of known allergy to some of the components do not use it.

Side effects

In deep cavities we suggest to use a liner in order to avoid pulpal reactions.

Materials to be avoided

Materials containing phenolics (like eugenol) and self etching primers could inhibit composite curing. Avoid the use of these materials as liners. This also applies to self-etching primers, which could prevent a proper curing of the cement. Therefore, the use of such primers should be foregone. We recommend the use of EnaBond.

DIRECTIONS FOR USE

1. Post cementation

A. Complete the endodontic therapy, and prepare the tooth for the restoration by evaluating if a post is necessary, based on residual quantity of intact tooth. The loss of one or both interproximal walls in anteriors, of one or both mesial or distal walls in posteriors, suggest you to use a post, in order to better stabilize the restoration.

B. Make a radiograph to determine the appropriate diameter and depth of post space preparation. The post diameter should be slightly higher than the endodontic preparation. The post should be placed at same radicular depth as the height of dentinal body (for composite aesthetic restoration) or of post and core. We recommend the use of rubber dam.

C. Remove gutta-percha to pre-planned depth with a Gates-Glidden

drill, Peeso reamer and/or hot instrument. Radiographic verification is recommended.

D. Select the post drill diameter (like EnaPost) which corresponds to the last Gates-Glidden drill to begin preparing the post space.

Use Drills with a slow-speed contra-angle with water spray in order to avoid tissues overheating. The drill must be kept in continuous clockwise rotation until it has been completely removed from the tooth. This will minimize the risk of the drill seizing in the post space. Remove any tooth debris from canal by irrigating post space with water spray. Sequentially step up to the next larger post drill until the pre-planned diameter and depth are achieved.

E. Select the post (like EnaPost) that corresponds to the last drill used to prepare the post space and insert post into post space. For a better fitting shorten the post as necessary from either the apical or occlusal end (as clinical judgement dictates) with a Diamond separating disc, wetting the post with water to avoid fibers overheating.

F. Remove the post and gently clean it with an alcohol wipe and apply on the surface a mixture of EnaBond with EnaBond Catalyst (one drop + one drop). A silane (EnaEtch silane) can also be used to increase the cement adhesion to the post.

G. Etch the cavity (coronal and radicular preparation) with EnaEtch 37% phosphoric acid for 2 minutes. For ideal adhesion, before etching, a microblasting on cavity surfaces is recommended in order to clean and eliminate endodontic material debris.

H. Wash accurately the canal with a syringe and water to remove completely the acid. Suck water and dry the canal with paper point; do not dry with air in order to maintain the dentine wet.

I. Apply a mixture of EnaBond and EnaBond Catalyst in the cavity and in the canal. The adhesive should be wiped on the surface with a disposable microbrush or with a paper point. **Attention:** Ensure that microbrush reaches into the depths of the canal and that the Bonder is evenly rubbed in everywhere. Microbrush should not touch the surface or possibly get jammed. Dry with air to eliminate water and solvent residuals. Insert the post to check the canal and better push the adhesive in dentine tubules.

J. Apply a disposable automixing tip on the syringe and press EnaCem HF dual composite cement out of the syringe: the two components will be automatically mixed. Always discard from the syringe the first ca. 0,5 g of material. Place the cement directly into the canal filling it from the bottom of the cavity to the surface using a disposable endodontic applicator (EnaCem HF orange tip); do not use any lentulo or rotating instrument. Apply some cement on the post surface and insert slowly the post to full depth, allowing excess cement to vent. Find the right position and, maintaining the post in position, wipe away any excess of cement. **Note:** after each use remove the mixing tip and apply again the cap to the syringe.

K. Light cure for 60 seconds and proceed to the restoration. Use a Mycro-Hybrid aesthetic composite (as Enamel Plus HFO), for aesthetic direct restoration, or Ena Cem for post and core. **Attention:** self curing in deep areas will be completed after 5-6 min. EnaCem HF working time is 3-4 min.

2. Core build-up

A. If you don't have already etched also the coronal surface as described in 1.G, etch the dentine surface with EnaEtch 37% phosphoric acid for 1 minute. Wash accurately the canal with a syringe to remove completely the acid. Suck water and dry the canal with paper point; do not dry with air in order to maintain the dentine wet.

B. Apply on the dentine surface the bonding (a dual bonding is not necessary in this case, but can be used if the core is over 2,5 mm thick: in this case follow the bonding instructions as from 1.I). The adhesive should be wiped on the surface with a disposable microbrush. Dry with air to eliminate water and solvent residuals. Light cure the adhesive; if you use EnaBond without catalyst light cure for 60 sec. and then repeat step B and light cure again.

C. Apply a disposable automixing tip and press EnaCem HF dual composite cement out of the syringe: the two components will be automatically mixed. Always discard from the syringe the first ca. 0.5 g of material. Place the cement directly on the dentine surface over the post. If you use a coreformer, fill it and put it on the post, finding the right position. Wipe away any excess cement. **Note:** after each use remove the mixing tip and apply again the cap to the syringe.

D. Light cure for 60 seconds. **Attention:** self curing in deep areas will be completed after 5-6 min. EnaCem HF working time is 3-4 min.

E. Finish and polish and take the impression.

3. Composite crowns & bridges, veneer, inlay and onlay luting

A. Remove the temporary appliance and clean the cavity. Try-in the appliance carefully and proceed with eventual corrections. You can use the special Try-in EnaCem HF bodies available in the same colours of EnaCem HF cement in order to find the right colour: the EnaCem HF Try-in will not cure. After finding the right colour wash away the try-in material.

B. Apply the rubber dam. Clean the surface of the preparation with alcohol and sandblast it.

C. Etch the tooth surface with EnaEtch 37% phosphoric acid for 1 minute.

D. Apply on the preparation the mixture of EnaBond and EnaBond Catalyst. The adhesive should be wiped on the surface with a disposable microbrush. Dry with air to eliminate water and solvent residuals.

E. Sandblast the internal part of the composite appliance and clean it with alcohol; apply EnaBond without curing it but dry it carefully.

F. Apply a small amount of EnaCem HF in the colour you have chosen in the internal side of the appliance to be luted, position it on the tooth and condense it mechanically or manually. Remove composite excess and light-cure for 60 seconds from each side of the tooth. **Attention:** self curing in deep areas will be completed after 5-6 min. EnaCem HF working time is 3-4 min.

G. Check the occlusion, finish and polish with Enamel plus Shiny system, using burs, strips and diamond pastes.

4. Ceramic & Metal crowns & bridges, veneer, inlay and onlay luting

A. Remove the temporary appliance and clean the cavity. Try-in the appliance carefully and proceed with eventual corrections. You can use the special Try-in EnaCem HF bodies available in the same colours of EnaCem HF cement in order to find the right colour: the EnaCem HF Try-in will not

- cure. After finding the right colour wash away the Try-in material with ethanol.
- B. Apply the rubber dam. Clean the surface of the preparation with alcohol and sandblast it.
- C. Etch the tooth surface with EnaEtch 37% phosphoric acid for 1 minute.
- D. Apply on the preparation the mixture of EnaBond and EnaBond Catalyst. The adhesive should be wiped on the surface with a disposable microbrush. Dry with air to eliminate water and solvent residuals.
- E. Always condition the contact surfaces of the restoration according to the manufacturer's instructions. In case of ceramic we suggest to etch the inner part of the application with hydrofluoric acid, wash it perfectly and then apply a silane. On metal can be used a special metal primer system like Tender Bonding that is ideal also for synthetic ceramic in combination with Tender paste opaque clear.
- F. Apply a small amount of EnaCem HF in the colour you have chosen in the internal side of the appliance to be luted, position it on the tooth and condense it mechanically or manually. Remove composite excess and cure for 60 seconds from each side of the tooth. **Attention:** self curing in deep areas will be completed after 5-6 min. EnaCem HF working time is 3-4 min.
- G. Check the occlusion, finish and polish.

5. Liner in posterior restoration

EnaCem HF can also be used in posterior restoration as a liner without light curing. In this way it will self cure creating a lower shrinkage towards the base of the tooth (instead towards the light as in light curing), without any stress on the tooth walls. Further the material will be harder and more resistant. No bubbles (typical for flowable composites) will be observed.

A. Etch and bond as usually.

B. Apply a small amount of EnaCem HF in the colour you have chosen directly on the preparation (original Enamel Plus HFO shade-guide) and let it self cure for 5-6 min. **Attention:** EnaCem HF working time is 3-4 min.

C. Proceed with the application of the light curing composite as usual, or finish the surface of EnaCem HF and take an impression in case an indirect technique has been chosen.

Curing information

For light curing it's necessary to use a light-curing unit with a spectrum of 350 - 500 nm. The required physical results can be reached only if using multi walls reflecting unit. For this reason we suggest a periodical check of the light intensity following the manufacturer's instructions.

Material can also self cure without light. Self curing in deep areas will be completed after 5-6 min. EnaCem HF working time is 3-4 min.

USE AND STORAGE

Store the product between 2°C and 10°C (between 36°F and 50°F) in a refrigerator. Do not use the product after the expiration date (see label on syringe). Due to hygienic reasons application tips should be used only once. The luting composite shows no chemical affinity to metals.

Use the material at room temperature. Medical device, for dental use only: keep away from children. After each use remove the mixing tip and apply again the cap to the syringe. Avoid direct sunlight exposure.

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