

AMERICAN ORTHODONTICS announces the addition of the Herbst™ Kit to its broad and innovative line of functional appliances. Available in one patient kits or as individual components, the Herbst is one of the most popular appliances used today to treat Class II corrections.

FEATURES

- Hex head screws for ease of use and maximum patient comfort; compatible with either straight or offset rod style Herbst.
- Offset rods provide a lower profile for special situations.
- Advancers to reactivate mandibular advancement are available in five lengths from 1 mm to 5 mm.
- Pivot casings for safe and secure solders every time. Two styles available: 'Regular' for soldering onto bands and crowns, and 'Flat' for soldering onto wire extensions.

HERBST
appliance

